

# MK-8000D



# MK-9000A/B/C

MK-9000B



## MATERNAL / FETAL MONITOR

### Features & Benefits:

- 15 inch high resolution color TFT LCD screen.
- Compact design, large screen for easy observation.
- Pace-maker detection, real time S-T segment analysis.
- Striking alarm light, 3-level audio/visual alarm.
- Multi-lead ECG wave forms display, capture dynamic wave forms.
- Large volume of tabular and graphic trends information storage and easy to recall.
- Efficient resistance to interference of defibrillator and HF knife.
- UP to 4 hours working capacity with built-in rechargeable lithium battery.
- Support solution: wall mounting, mobile trolley.
- Wire/wireless central monitoring system supported.

MK-8000E



10.4 inch color TFT LCD screen

MK-8000F



8 inch color TFT LCD display

#### Standard Configuration

SpO2, PR, NIBP, ECG, RESP, TEMP

#### Optional Configuration

Nellcor SpO2, Dual IBP, ETCO2, Touch screen, Thermal printer

### Introduction:

- \* 12.1" high resolution color TFT LCD display. \* Simple and friendly operating display interface..
- \* 9 Standard parameters: ECG, RESP, NIBP, SPO2, 1-TEMP, PR/HR, FHR, TOCO, FM.

### Features & Benefits:

- 9 Standard parameters: ECG, RESP, NIBP, SPO2, 1-TEMP, PR/HR, FHR, TOCO, FM.
- 12.1"high resolution color TFT LCD display. Simple and friendly operating display interface.
- 9-wafer impulse ultrasound transducer, configured with special HR self correlation reading algorithm to
  - ensure correct and reliable fetal heart rate (FHR) measurement
- Large built-in chargeable battery provides convenience for patients' moving.
- Playback and browse function for history waveform and monitor data record.
- Optional printing output function, alarm triggers printing.
- Auto double alarm with audible and visible signals.
- Anti-defibrillator, anti-interference from high frequency electric knife.
- Pace-maker detection.
- A full-synchronistic leads, multi-channel ECG display.
- Indication of incorrect manipulation and malfunction analysis.
- Unique accessory management cabinet.

MK-9000C



MK-9000A



Model	FHR	TOCO	FM
MK-9000A	√	√	√
MK-9000C	√	√	√